

SCISSORS



Kazan Medical Instruments Plant



Tips 0,8 mm wide
Blades 16 mm long

S-0165

**WESTCOTT SCISSORS
CURVED, STAINLESS STEEL**

- Overall length 135 mm; 5,3 inches
- Flat handle 8 mm wide
- Mat finish



Tips 0,25 mm wide
Blades 8,5 mm long

S-0160

**CASTROVIEJO MINI CORNEAL
SCISSORS CURVED,
STAINLESS STEEL**

- Overall length 125 mm; 4,9 inches
- Round handle 7 mm dia
- Mat finish

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

SCISSORS

RS



Tips 0,8 mm wide
Blades 16 mm long

S-0000



Tips 0,4 mm wide
Blades 16 mm long

S-0000.01



WESTCOTT SCISSORS

CURVED, STAINLESS STEEL

Overall length 120 mm; 4,7 inches

Flat handle 8 mm wide

Mat finish



Tips 0,6 mm wide
Blades 14 mm long

S-0140



Tips 0,6 mm wide
Blades 14 mm long

S-0140.01



RIGHT

LEFT

WESTCOTT SCISSORS CURVED, STAINLESS STEEL

Overall length 120 mm; 4,7 inches

Flat handle 8 mm wide

Mat finish

STAINLESS

SCISSORS



Kazan Medical Instruments Plant



Tips 0,3 mm wide
Blades 8 mm long

S-0026**STRAIGHT**

Tips 0,3 mm wide
Blades 9,5 mm long

S-0027**CURVED**

Tips 0,3 mm wide
Blades 8 mm long
Angled 15 degrees

S-0029**ANGLED TO FLAT**
**VANNAS SCISSORS,
STAINLESS STEEL**

- Overall length 118 mm; 4,6 inches
- Flat handle 7 mm wide
- Mat finish



Tips 0,15 mm wide
Blades 6 mm long

S-0130**STRAIGHT**

Tips 0,15 mm wide
Blades 6 mm long

S-0130.01**CURVED**

Tips 0,2 mm wide
Blades 5,5 mm long
Angled 20 degrees

S-0143**ANGLED TO FLAT**
**VANNAS TYPE SCISSORS,
STAINLESS STEEL**

- Overall length 105 mm; 4,2 inches
- Round handle 8 mm dia with two hinge joints
- Mat finish

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

SCISSORS SCISS

ORS



Tips 0,3 mm wide
Blades 8,5 mm long

S-0141

● **CASTROVIEJO SCISSORS**
● **CURVED,**
● **STAINLESS STEEL**

Overall length 110 mm; 4,3 inches
Flat handle 8 mm wide
Mat finish



Tips 0,65 mm wide
Blades 14 mm long

S-0064

● **CASTROVIEJO UNIVERSAL**
● **SCISSORS CURVED,**
● **STAINLESS STEEL**

Overall length 105 mm; 4,1 inches
Flat handle 7 mm wide
Mat finish



Tips 0,4 mm wide
Blades 14 mm long

S-0136



Tips 0,8 mm wide
Blades 14 mm long

S-0137

● **WESTCOTT SCISSORS CURVED,**
● **STAINLESS STEEL**

Overall length 110 mm; 4,3 inches
Flat handle 7 mm wide
Mat finish

STAINLESS

SCISSORS



Kazan Medical Instruments Plant



Tips 0,25 mm wide
Blades 10 mm long
Angled 45 degrees

S-0074

ANGLED TO FLAT



Tips 0,3 mm wide
Blades 8,5 mm long

S-0075

CURVED



Tips 0,3 mm wide
Blades 10 mm long
Angled 40 degrees

S-0122

ANGLED TO SIDE

SCISSORS, STAINLESS STEEL

- Overall length 90 mm; 3,5 inches
- Flat handle 6 mm wide
- Mat finish



Tips 0,3 mm wide
Blades 11,5 mm long

S-0066

CURVED



Tips 0,65 mm wide
Blades 12,5 mm long

S-0079

STRAIGHT

SCISSORS FOR CONJUNCTIVA, STAINLESS STEEL

- Overall length 100 mm; 3,9 inches
- Flat handle 7 mm wide
- Mat finish

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

RS



Tips 0,2 mm wide
Blades 12 mm long

S-0040



SHARP
Tips 0,2 mm and
0,5 mm wide
Blades 12 mm long

S-0040.01



COMBINED
Tips 0,5 mm
Blades 12 mm long

S-0040.02

BLUNT



**NOYES SCISSORS STRAIGHT,
STAINLESS STEEL**

Overall length 106 mm; 4,1 inches
Flat handle 7 mm wide
Mat finish



Tips 0,2 mm wide
Blades 8 mm long

S-0043



SHARP
Tips 0,2 mm and
0,5 mm wide
Blades 8 mm long

S-0043.01



COMBINED
Tips 0,5 mm wide
Blades 8 mm long

S-0043.02

BLUNT



**NOYES SCISSORS CURVED,
STAINLESS STEEL**

Overall length 100 mm; 3,9 inches
Flat handle 7 mm wide
Mat finish

SCISSORS SCISSO

STAINLESS

SCISSORS SCISSORS



Kazan Medical Instruments Plant



Tips 0,5 mm wide
Blades 10 mm long
Angled 50 degrees

S-0072



Tips 0,5 mm wide
Blades 10 mm long
Angled 50 degrees

S-0073



RIGHT



LEFT

**TROUTMAN - CASTROVIEJO
CORNEAL SECTION SCISSORS,
STAINLESS STEEL**

- Overall length 100 mm; 3,9 inches
- Flat handle 6 mm wide
- Mat finish



Tips 0,5 mm wide
Blades 13 mm long
Angled 50 degrees

S-0134



Tips 0,5 mm wide
Blades 13 mm long
Angled 50 degrees

S-0134.01



RIGHT



LEFT

**CASTROVIEJO CORNEAL
SECTION SCISSORS,
STAINLESS STEEL**

- Overall length 105 mm; 4,1 inches
- Flat handle 7 mm wide
- Mat finish

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

SCISSORS

RS



Tips 0,3 mm wide
Blades 7 mm long
Angled 45 degrees
S-0069



Tips 0,3 mm wide
Blades 7 mm long
Angled 45 degrees
S-0069.01



RIGHT



LEFT

- **KATZIN CORNEAL TRANSPLANT**
- **SCISSORS, STAINLESS STEEL**
- Overall length 100 mm; 3,9 inches
- Round handle 8 mm dia
- Mat finish



Tips 0,4 mm wide
Blades 9 mm long
Angled 30 degrees
S-0080



Tips 0,4 mm wide
Blades 9 mm long
Angled 30 degrees
S-0081



RIGHT



LEFT

- **CASTROVIEJO CORNEAL SECTION**
- **SCISSORS, STAINLESS STEEL**
- Overall length 100 mm; 3,9 inches
- Flat handle 6 mm wide
- Mat finish

STAINLESS

SCISSORS



Kazan Medical Instruments Plant



Tips 0,5 mm wide
Blades 12 mm long
Angled 45 degrees

S-0164

**ANTERIOR CAPSULOTOMY
SCISSORS ANGLED,
STAINLESS STEEL**

- Overall length 95 mm; 3,7 inches
- Flat handle 7 mm wide
- Mat finish



Tips 0,15 mm wide
Blades 6 mm long

S-0135

Tips 0,15 mm wide
Blades 6 mm long

S-0135.01

Tips 0,15 mm wide
Blades 6 mm long

S-0146

Tips 0,15 mm wide
Blades 6 mm long

S-0146.01

Tips 0,3 mm wide
Blades 7 mm long

S-0139

**WESTCOTT SCISSORS
CURVED, STAINLESS STEEL**

- Overall length 84 mm; 3,3 inches
- Flat handle 7 mm wide
- Mat finish

**CAPSULOTOMY SCISSORS
SHARP, STAINLESS STEEL**

- Overall length 82 mm; 3,2 inches
- Flat handle 8 mm wide
- Mat finish

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

RS



VANNAS SCISSORS,
STAINLESS STEEL
Flat handle 7 mm wide
Mat finish



Tips 0,4 mm wide
Blades 9,5 mm long
S-0044

Overall length 84 mm; 3,3 inches



Tips 0,4 mm wide
Blades 9,5 mm long
S-0045

Overall length 84 mm; 3,3 inches



Tips 0,4 mm wide
Blades 9,5 mm long
Angled 45 degrees
S-0046

Overall length 84 mm; 3,3 inches




Tips 0,4 mm wide
Blades 9,5 mm long
Angled 50 degrees
S-0047

Overall length 84 mm; 3,3 inches



Tips 0,2 mm wide
Blades 5,5 mm long
S-0044.01

Overall length 80 mm; 3,1 inches



Tips 0,2 mm wide
Blades 5,5 mm long
S-0045.01

Overall length 80 mm; 3,1 inches



Tips 0,15 mm wide
Blades 6 mm long
S-0046.01

Overall length 80 mm; 3,1 inches



Tips 0,2 mm wide
Blades 5,5 mm long
Angled 50 degrees
S-0047.01

Overall length 80 mm; 3,1 inches

SCISSORS

SCISSO

STAINLESS

SCISSORS



Kazan Medical Instruments Plant



Tips 0,3 mm wide
Blades 9 mm long

S-0018

MICRO

Overall length 100 mm; 3,9 inches



Tips 0,3 mm wide
Blades 12 mm long

S-0065

Overall length 105 mm; 4,1 inches



Tips 0,6 mm wide
Blades 13,5 mm long

S-0070

UNIVERSAL

Overall length 108 mm; 4,2 inches

**CASTROVIEJO SCISSORS
CURVED, STAINLESS STEEL**

- Round handle 8 mm dia
- Mat finish



Tips 0,3 mm wide
Blades 10 mm long
Angled 45 degrees

S-0067

ANGLED TO FLAT

Overall length 82 mm; 3,2 inches



Tips 0,3 mm wide
Blades 8,5 mm long

S-0068

CURVED

Overall length 84 mm; 3,3 inches

SCISSORS, STAINLESS STEEL

- Round handle 6 mm dia
- Mat finish

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

RS



Tips 0,2 mm wide
Blades 6 mm long

S-0163

Overall length 112 mm;
4,4 inches

Tips 0,2 mm wide
Blades 7 mm long
Angled 30 degree

S-0163.01

Overall length 113 mm;
4,4 inches

Tips 0,2 mm wide
Blades 8 mm long

S-0163.02

Overall length 114 mm;
4,5 inches

Tips 0,2 mm wide
Blades 6 mm long

S-0163.03

SCISSORS WITH TITANIUM HANDLE

Flat handle 7,5 mm wide
with two hinge joints
Mat finish



Tips 0,4 mm wide
Blades 9,0 mm long
Jaws 1,0 mm long
Angled 45 degree

S-0171



SCISSORS-FORCEPS WITH TITANIUM HANDLE

Overall length 113 mm; 4,4 inches
Flat handle 7,5 mm wide with two hinge joints
Mat finish

Overall length 112 mm;
4,4 inches

SCISSORS

SCISSO

STAINLESS

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant



SPECULUM


TITANIUM



SPECULUM SPECULUM



Kazan Medical Instruments Plant




R-0000
with movable supports




R-0002
without movable supports

EYE SPECULUM, TITANIUM

- Overall length 80 mm; 3,1 inches
- Close 15 mm wire blades



R-0001
RIGHT



R-0001.01
LEFT

MURDOCK EYE SPECULUM, TITANIUM

- Overall length 70 mm; 2,7 inches
- Close 14 mm wire blades



Kazan Medical Instruments Plant

SPECULUM SPECCULUM



**BARRAQUER EYE LID
SPECULUM**

**RT-0001
ADULT**

- Overall length 43 mm; 1,7 inches
- Close 15 mm wire blades
- Dull finish

**RT-0001.01
PAEDIATRIC**

- Overall length 40 mm; 1,6 inches
- Close 10 mm wire blades
- Dull finish




**KRATZ-BARRAQUER EYE LID
SPECULUM**

**RT-0002.01
ADULT**

- Overall length 44 mm; 1,7 inches
- Open 14 mm wire blades
- Dull finish

**RT-0002
PAEDIATRIC**

- Overall length 38 mm; 1,5 inches
- Open 10 mm wire blades
- Dull finish



**BARRAQUER EYE LID
SPECULUM HEAVY ADULT**

**RT-0003
ADULT**

- Overall length 50 mm; 2,0 inches
- Close 14 mm wire blades
- Dull finish



**KRATZ-BARRAQUER EYE LID
SPECULUM**

**R-0003
ADULT**

- Overall length 41 mm; 1,6 inches
- Close 14 mm wire blades
- Dull finish

SPECULUM

SPECULUM



Kazan Medical Instruments Plant



R-0010
ADULT

- Overall length 80 mm; 3,1 inches
- Open 16 mm wire blades
- Angled 45 degrees

R-0010.01
PAEDIATRIC

- Overall length 77 mm; 3,0 inches
- Open 12 mm wire blades
- Angled 45 degrees

R-0010.04
PAEDIATRIC

- Overall length 60 mm; 2,3 inches
- Open 10 mm wire blades
- Angled 45 degrees



R-0010.02
ADULT

- Overall length 70 mm; 2,7 inches
- Open 16 mm wire blades
- Angled 45 degrees

R-0010.03
PAEDIATRIC

- Overall length 68 mm; 2,7 inches
- Open 12 mm wire blades
- Angled 45 degrees

R-0010.05
PAEDIATRIC

- Overall length 64 mm; 2,5 inches
- Open 10 mm wire blades
- Angled 45 degrees



R-0012
ADULT

- Overall length 88 mm; 3,4 inches
- Open 16 mm wire blades
- Angled 45 degrees

R-0012.01
PAEDIATRIC

- Overall length 79 mm; 3,1 inches
- Open 12 mm wire blades
- Angled 45 degrees

**SPECULUM ADJUSTABLE ANGLED,
TITANIUM**



Kazan Medical Instruments Plant

SPECULUM SPECULUM




**R-0011
ADULT**

- Overall length 67 mm; 2,7 inches
- 15 mm solid blades
- Angled 45 degrees

**R-0011.01
PAEDIATRIC**

- Overall length 65 mm; 2,5 inches
- 11 mm solid blades
- Angled 45 degrees



**R-0013.02
ADULT**

- Overall length 69 mm; 2,7 inches
- 15 mm solid blades
- Angled 45 degrees

**R-0013.03
PAEDIATRIC**

- Overall length 67 mm; 2,7 inches
- 11 mm solid blades
- Angled 45 degrees



**R-0013
ADULT**

- Overall length 82 mm; 3,2 inches
- 15 mm solid blades
- Angled 45 degrees

**R-0013.01
PAEDIATRIC**

- Overall length 80 mm; 3,1 inches
- 11 mm solid blades
- Angled 45 degrees

EYE SPECULUM ADJUSTABLE ANGLED,
TITANIUM

SPECULUM SPECULUM



Kazan Medical Instruments Plant



SMP-0015

**SCHAPIRA SPECULUM,
TITANIUM**


- Overall length 94 mm; 3,7 inches
- Open 14 mm wire blades
- Angled 40 degrees



R-0032

EYE SPECULUM, TITANIUM

- Overall length 80 mm; 3,1 inches
- Open 15 mm wire blades
- Angled 40 degrees



R-0033

**McBRATHEY ASPIRATING
SPECULUM, STAINLESS STEEL**

- Overall length 42 mm; 1,7 inches
- Open 14 mm wire blades
- Angled 30 degrees
- Dull finish



Kazan Medical Instruments Plant



**R-0027
SMALL**

- Close 14 mm wire blades

**R-0027.01
MEDIUM**

- Close 16 mm wire blades


**R-0027.02
LARGE**

- Close 18 mm wire blades

**CASTROVIEJO EYE
SPECULUM, STAINLESS STEEL**

- Overall length 94 mm; 3,7 inches

SPECULUM SPECULUM



**SPRING EYE
SPECULUM ADULT,
TITANIUM**

R-0005

- Overall length 60 mm; 2,4 inches
- Open 15 mm wire blades



**BARRAQUER EYE LID
SPECULUM,
STAINLESS STEEL**

**R-0008
PAEDIATRIC**

- Overall length 41 mm; 1,6 inches
- 11 mm solid blades
- Dull finish

**R-0008.01
ADULT**

- Overall length 42 mm; 1,65 inches
- 15 mm solid blades
- Dull finish

SPECULUM SPECULUM



Kazan Medical Instruments Plant



R-0031

EYE SPECULUM WITH ASPIRATION

- Overall length 77 mm; 3,0 inches
- Open 12 mm wire stainless steel blades
- Angled 45 degrees
- Titanium mechanism



EYE SPECULUM WITH ASPIRATION

RA-100
ADULT

- Overall length 88 mm; 3,4 inches
- Open 16 mm wire stainless steel blades
- Angled 45 degrees
- Titanium mechanism

RA-100.01
PEDIATRIC

- Overall length 79 mm; 3,1 inches
- Open 12 mm wire stainless steel blades
- Angled 45 degrees
- Titanium mechanism



OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant



ROD
INSTRUMENTS with

TITANIUM
HANDLES



ROD INSTRUMENTS



Kazan Medical Instruments Plant



Button 0,45 mm dia

E-0000

**GRAETHER COLLAR
BUTTON HOOK**



Tip 0,4 mm wide
0,2 mm dia

E-0000.06

**SINSKEY
IRIS and IOL HOOK**



Tip 0,6 mm wide
0,2 mm dia

E-0000.09

BONN IRIS HOOK



Tip 0,65 mm wide
0,65 mm long

E-0000.07

**HIRSHMAN IRIS
HOOK and LENS
MANIPULATOR**

STRAIGHT MICROHOOKS

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Button 0,45 mm dia
Angled 30 degrees

E-0000.12

**GRAETHER COLLAR
BUTTON HOOK**



Tip 0,4 mm wide
0,2 mm dia
Angled 30 degrees

E-0000.01

**SINSKEY
IRIS and IOL HOOK**



Tip 0,6 mm wide
0,2 mm dia
Angled 30 degrees

E-0000.03

BONN IRIS HOOK



Tip 0,65 mm wide
0,65 mm long
Angled 30 degrees

E-0000.02

**HIRSHMAN IRIS
HOOK and LENS
MANIPULATOR**

ANGLED MICROHOOKS

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Kazan Medical Instruments Plant

ROD INSTRUMENTS ROD INSTRUMENTS



Tip 0,85 mm wide

E-0000.13

OSHER & Y[HOOK



Tip 0,15 mm dia

E-0071

**FENZL INSERTION
HOOK**



Hook 1,0 mm wide

E-0013

**SHEPARD IRIS
HOOK**



Tip 0,85 mm wide
0,15 mm dia

E-0000.04

**RENTSCH LENS
MANIPULATOR**



Tip 0,85 mm wide
Angled 30 degrees

E-0000.05

OSHER & Y[HOOK



Tip 0,15 mm dia
Angled 40 degrees

E-0071.01

**FENZL INSERTION
HOOK**



Tip 0,5 mm wide
0,15 mm dia
Angled 30 degrees

E-0001

MICROHOOK RIGHT MICROHOOK LEFT



Tip 0,5 mm wide
0,15 mm dia
Angled 30 degrees

E-0001.01

STRAIGHT MICROHOOKS

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips

ANGLED MICROHOOKS

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips

ROD INSTRUMENTS



Kazan Medical Instruments Plant



Tip 1,5 mm wide

E-0012**RIGHT**

Tip 1,5 mm wide

E-0012.01**LEFT****TWIST HOOK SMALL**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips

Tip 5,0 mm wide
0,6 mm dia**E-0003****MAUMENEE IRIS HOOK**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tip

Tip 0,35 mm dia
Angled 45 degrees**E-0005****IOL MANIPULATION HOOK**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tip

Sharp tip
2,0 mm wide**E-0002****LARGE**Sharp tip
1,5 mm wide**E-0002.01****SMALL****DOUBLE FIXATION HOOK,
STRAIGHT**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant



MICROKNIVES, ANGLED

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Tip 0,9 mm dia
Blade 1,5 mm long
Angled 45 degrees

K-0001

LENS CRUSHING KNIFE



Tip 0,3 mm wide
Angled 60 degrees

K-0002

CYSTOTOME



Tip 1,0 mm wide
Blade 1,8 mm long
Angled 45 degrees

K-0003

DRISDALE LENS CRUSHING KNIFE



Tip 1,0 mm wide
Blade 1,8 mm long
Angled 45 degrees

K-0005

DRISDALE LENS CRUSHING KNIFE CURVED

ROD INSTRUMENTS

ROD INSTRUMENTS



Blade 1,5 mm long

K-0010.04



Blade 1,3 mm long

K-0010.05



Blade 1,0 mm long

K-0010.06

LENS CRUSHING KNIFE ANGLED, TITANIUM

- Overall length 120 mm; 4,7 inches
- Round handle 5,0 mm dia
- Angled 45 degrees
- Stainless steel tips

TITANIUM
ROD INSTRUMENTS

ROD INSTRUMENTS

Kazan Medical Instruments Plant



Tip 1,5 mm dia

K-0031

Tip 1,75 mm dia

K-0031.01

Tip 2,0 mm dia

K-0031.02

Tip 2,5 mm dia

K-0031.03

Tip 3,5 mm dia

K-0031.04**MEYERHOEFFER CHALAZION
CURETTES**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Kazan Medical Instruments Plant



Blade 1,0 mm long

K-0000

SATO KNIFE



Tip 1,8 mm wide

K-0017

**FOREIGN BODY
REMOVAL SPUD**

Lance shaped gently
curved tip



Tip 0,6 mm wide

K-0016

GULF CLUB SPUD

Small foreign body
angled tip with edges



Tip 1,2 mm wide

K-0018

**DAVIS FOREIGN BODY
REMOVAL SPUD**

Straight sharp round blade



Blade 1,2 mm long

K-0019

FRANCIS SPUD

Small foreign body removal flag
shaped angled with sharp edges

STRAIGHT MICROSPUDS

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips

ROD INSTRUMENTS



Kazan Medical Instruments Plant



Tip 4,0 mm wide
Loop 14 mm long

L-0002**WILDER LOOP**

Tip 6,0 mm wide
Angled 30 degrees

L-0000**WEBER LOOP****LENS LOOPS**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips

**M-0002**

- Blades external diameter 12 mm

- M-0002 3 blades
- M-0002.01 4 blades
- M-0002.02 5 blades
- M-0002.05 7 blades

**M-0003**

- Blades external diameter 14 mm

- M-0003 4 blades
- M-0003.01 6 blades
- M-0003.02 8 blades
- M-0003.03 10 blades
- M-0009.04 12 blades

**THORNTON RADIAL MARKERS
LOW PROFILE**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel head



Kazan Medical Instruments Plant



M-0007

**BORES AXIS MARKER
USED WITH MENDEZ
DEGREE GAUGE**

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel head



Z-0008

**MENDEZ DEGREE GAUGE
FOR MAKING AXIS OF
ASTIGMATISM, TITANIUM**

- Overall length 135 mm; 5,3 inches
- Round handle 5,0 mm dia
- 12 mm internal dia
- Angled 45 degrees

ROD INSTRUMENTS



Kazan Medical Instruments Plant



M-0008
OUTSIDE OPTIC
ZONE MARKER



M-0009
INSIDE OPTIC
ZONE MARKER,
LOW PROFILE



M-0010
THORTOH OPTIC
ZONE MARKER
WITH CENTERING PIN,
LOW PROFILE



M-0008	M-0008.05	M-0008.10	M-0008.15
M-0009 2,50	M-0009.05 3,75	M-0009.10 5,00	M-0009.15 6,5
mm	mm	mm	0 mm
M-0008.01	M-0008.06	M-0008.11	M-0008.16
M-0009.01 2,7	M-0009.06 4,00	M-0009.11 5,25	M-0009.16 7,0
5 mm	mm	mm	0 mm
M-0008.02	M-0008.07	M-0008.12	M-0008.17
M-0009.02 3,0	M-0009.07 4,25	M-0009.12 5,50	M-0009.17 7,5
0 mm	mm	mm	0 mm
M-0008.03	M-0008.08	M-0008.13	M-0008.18
M-0009.03 3,2	M-0009.08 4,50	M-0009.13 5,75	M-0009.18 8,0
5 mm	mm	mm	0 mm
M-0008.04	M-0008.09	M-0008.14	M-0008.19
M-0009.04 3,5	M-0009.09 4,75	M-0009.14 6,00	M-0009.19 8,5
0 mm	mm	mm	0 mm
M-0010.04	M-0010.09	M-0010.14	M-0010.19

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Kazan Medical Instruments Plant

ROD INSTRUMENTS



R-0006

- Blade 10 mm wide
- Overall length 120 mm; 4,7 inches



R-0006.01

- Blade 11 mm wide
- Overall length 120 mm; 4,7 inches



R-0006.02

- Blade 13 mm wide
- Overall length 125 mm; 4,9 inches



R-0006.03

- Blade 15 mm wide
- Overall length 126 mm; 4,9 inches



R-0006.04

- Blade 17 mm wide
- Overall length 126 mm; 4,9 inches

DESMARRES LID RETRACTORS

- Round titanium handle 5,0 mm dia
- Stainless steel tips

ROD INSTRUMENTS



Kazan Medical Instruments Plant



Tips 3,5 mm wide
Angled 15 degree



Tip 2,2 mm dia
Angled 25 degrees

R-0007

IRIS RETRACTOR

- Overall length 150 mm; 5,9 inches
- Round titanium handle 6,0 mm dia
- Stainless steel tips



Tip 1,0 mm wide
0,3 mm dia

T-0014

CONVEX



Tip 1,0 mm wide
0,3 mm dia

T-0014.01

CONCAVE

BANANA-SHAPED EYE SPATULA

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5 mm dia
- Angled 30 degrees
- Stainless steel tips



Kazan Medical Instruments Plant

ROD INSTRUMENTS ROD INSTRUMENTS



Tip 0,5 mm wide
14 mm long
0,25 mm thickness
T-0000



Tip 1,0 mm wide
14 mm long
0,25 mm thickness
T-0000.01



Tip 1,5 mm wide
14 mm long
0,25 mm thickness
T-0000.02



Tip 0,8 mm wide
15 mm long
0,25 mm thickness
T-0001

**IOL NOTCHED
SPATULA**



Tip 0,25 mm wide
14 mm long
0,2 mm thickness
T-0006

**BARRAQUER
CYCLODIALYSIS
SPATULA**



Tip 0,5 mm wide
14 mm long
0,2 mm thickness
T-0006.01

**BARRAQUER
CYCLODIALYSIS
SPATULA**



Tip 0,75 mm wide
14 mm long
0,2 mm thickness
T-0006.02

**CASTROVIEJO
CYCLODIALYSIS
SPATULA**



Tip 1,0 mm wide
14 mm long
0,2 mm thickness
T-0006.03

**CASTROVIEJO
CYCLODIALYSIS
SPATULA**

SPATULAE, ANGLED

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips
- Angled 30 degrees

ROD INSTRUMENTS



Tip 4,0 mm wide

T-0003

KIMURA PLATINUM SPATULA

- Overall length 150 mm; 5,9 inches
- Round titanium handle 6,0 mm dia
- Platinum tip



Tip 0,5 mm wide
14 mm long
Angled 30 degrees



Tip 1,0 mm wide
14 mm long
Angled 30 degrees

T-0002

DOUBLE-ENDED EYE SPATULA

- Overall length 130 mm; 5,1 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Tip 0,25 mm wide
14 mm long
Angled 30 degrees



Tip 0,5 mm wide
14 mm long
Angled 30 degrees

T-0007

CASTROVIEJO DOUBLE-ENDED CYCLODIALYSIS SPATULA

- Overall length 120 mm; 4,7 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



Kazan Medical Instruments Plant

ROD INSTRUMENTS ROD INSTRUMENTS



Z-0000

- Overall length 135 mm;
5,31 inches
- 12 teeth



Z-0000.01

- 3/4 GLOBE
- Overall length 130 mm;
5,1 inches
- 9 teeth



Z-0000.02

- Overall length 135 mm;
5,3 inches
- 12 teeth



Z-0000.03

- Overall length 130 mm;
5,1 inches
- 9 teeth

THORNTON REVOLVING GLOBE FIXATION RINGS, TITANIUM

- Internal dia 16 mm (also available internal dia 13 mm)
- Round handle 5,0 mm dia
- To stabilize and rotate eyeglobe
- 3/4 globe can be used for AK,RK and clear corneal incision

ROD INSTRUMENTS



Kazan Medical Instruments Plant



Z-0006



Z-0006
HOFFER OPTIC
ZONE MARKER
WITH CENTERING
PIN, LOW PROFILE

Z-0006	2,00 mm	Z-0006.05	3,25 mm
Z-0006.01	2,25 mm	Z-0006.06	3,50 mm
Z-0006.02	2,50 mm	Z-0006.07	3,75 mm
Z-0006.03	2,75 mm	Z-0006.08	4,00 mm
Z-0006.04	3,00 mm	Z-0006.09	4,25 mm

Z-0006.10	4,50 mm	Z-0006.15	5,75 mm
Z-0006.11	4,75 mm	Z-0006.16	6,00 mm
Z-0006.12	5,00 mm	Z-0006.17	6,50 mm
Z-0006.13	5,25 mm	Z-0006.18	7,00 mm
Z-0006.14	5,50 mm	Z-0006.19	7,50 mm
		Z-0006.20	8,00 mm

- Overall length 120 mm; 4,72 inches
- Round titanium handle 5,0 mm dia
- Stainless steel tips



W-0000

EYE THERMOCAUTER

- Overall length 150 mm; 5,9 inches
- Tip 0,5 mm dia, oval working portion 5 mm dia
- Angled 45 degrees
- Round titanium handle 6,0 mm dia
- Stainless steel tip



Kazan Medical Instruments Plant

ROD INSTRUMENTS ROD INSTRUMENTS



Z-0000.07

- Overall length 120 mm; 4,7 inches
- 12 teeth



Z-0000.08

- 3/4 Globe
- Overall length 120 mm; 4,7 inches
- 9 teeth



Z-0000.09

- 1/2 Globe
- Overall length 115 mm; 4,5 inches
- 7 teeth

THORNTON CLOBE FIXATION RINGS, TITANIUM

- Round handle 5,0 mm dia
- Angled 45 degrees
- To stabilize and rotate eyeglobe
- Z-0000.08 and Z-0000.09 can be used for AK,RK and clear corneal incision



Z-0007

REVOLVING

- Overall length 130 mm; 5,1 inches



Z-0005

Z-0005

FIXATION

- Overall length 115 mm; 4,5 inches
- Angled 45 degrees

3/4 FIXATION RINGS & FIXATOR[, TITANIUM

- Round handle 5,0 mm dia
- Concentric atraumatic ring to stabilize and rotate eyeglobe
- Super instrument for AK,RK,clear corneal incision and trabeculoectomy

TITANIUM

ROD INSTRUMENTS

ROD INSTRUMENTS

Kazan Medical Instruments Plant



Tip 2,5 mm dia
Angled 45 degrees

D-0000

WILDER SCLERAL DEPRESSOR

- Overall length 125 mm; 4,9 inches
- Round titanium handle 6,0 mm dia
- Stainless steel tip



Tip 2,5 mm dia
1,0 mm wide



Tip 6 mm wide
2,5 mm dia


D-0001

SCHOCKET SCLERAL DEPRESSOR DOUBLE-ENDED

- Overall length 135 mm; 5,3 inches
- Round titanium handle 6,0 mm dia
- Stainless steel ti




Kazan Medical Instruments Plant



**B-0000
STRAIGHT**

- 13 mm conical jaws



**B-0008
ANGLED**

- 11 mm conical jaws
- Angled 45 degrees

CASTROVIEJO BLADE HOLDERS, TITANIUM

- Overall length 120 mm; 4,7 inches
- Round handle 6,0 mm dia
- TC jaws



B-0001

- 10 mm conical jaws

CASTROVIEJO BLADE HOLDER STRAIGHT, TITANIUM

- Overall length 95 mm; 3,7 inches
- Round handle 6,0 mm dia
- TC jaws

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

VARIOUS INSTRUMENTS VARIOUS INSTRUMENTS



M-0001

- Overall length 86 mm; 3,35 inches
- Measures up to 20 mm on either side



M-0001.01

- Overall length 76 mm; 3,00 inches
- Measures up to 15 mm on either side

CASTROVIEJO CALIPERS, TITANIUM

V-077

**TITANIUM INTRAOCULAR
INSTRUMENT STERILIZING BOX**



Removable cover with lock

- Dimensions 190 mm (7,48") x 130 mm (5,12") x 32 mm (1,26")




Kazan Medical Instruments Plant

K-0009

DONOR CORNEA PUNCH

- Supplied with set of trepans K-0008-0008.15
- Dimensions 100x80x120 mm



Z-0001

DONOR CORNEA HOLDER

- Dimensions 60x55 mm





Kazan Medical Instruments Plant



K-0006

EYE TREPANS CONIC

- Overall length 40 mm; 1,57 inches

K-0006	5,0 mm	K-0006.04	6,75 mm	K-0006.08	7,75 mm	K-0006.12	8,75 mm
K-0006.01	5,50 mm	K-0006.05	7,00 mm	K-0006.09	8,00 mm	K-0006.13	9,00 mm
K-0006.02	6,00 mm	K-0006.06	7,25 mm	K-0006.10	8,25 mm	K-0006.14	9,50 mm
K-0006.03	6,50 mm	K-0006.07	7,50 mm	K-0006.11	8,50 mm	K-0006.15	10,00 mm



K-0007

EYE TREPANS WITH TITANIUM HANDLE

- Overall length 66 mm; 2,6 inches
- Round handle 6 mm dia

K-0007	5,0 mm	K-0007.04	6,75 mm	K-0007.08	7,75 mm	K-0007.12	8,75 mm
K-0007.01	5,50 mm	K-0007.05	7,00 mm	K-0007.09	8,00 mm	K-0007.13	9,00 mm
K-0007.02	6,00 mm	K-0007.06	7,25 mm	K-0007.10	8,25 mm	K-0007.14	9,50 mm
K-0007.03	6,50 mm	K-0007.07	7,50 mm	K-0007.11	8,50 mm	K-0007.15	10,00 mm



K-0008

EYE TREPANS WITH TITANIUM HANDLE

- Overall length 55 mm; 2,2 inches
- Round handle 6 mm dia

K-0008	5,0 mm	K-0008.04	6,75 mm	K-0008.08	7,75 mm	K-0008.12	8,75 mm
K-0008.01	5,50 mm	K-0008.05	7,00 mm	K-0008.09	8,00 mm	K-0008.13	9,00 mm
K-0008.02	6,00 mm	K-0008.06	7,25 mm	K-0008.10	8,25 mm	K-0008.14	9,50 mm
K-0008.03	6,50 mm	K-0008.07	7,50 mm	K-0008.11	8,50 mm	K-0008.15	10,00 mm

OPHTHALMIC SURGICAL SET

S



Kozan Medical Instruments Plant

OPHTHALMIC
SURGICAL SETS

AO-1		
UNIVERSAL OPHTHALMIC SET		
1	S-0045	Vannas scissors, curved
2	S-0046	Vannas scissors angled to flat
3	S-0066	Scissors for conjunctiva curved
4	S-0043	Noyes scissors curved (sharp tips)
5	F-0224	Troutman-Barraquer corneal forcep colibri type, titanium
6	F-0225	Troutman-Barraquer corneal forcep colibri type, titanium
7	F-0257	Bonn iris forceps straight, titanium
8	F-0008	Catalano tying forceps curved, titanium
9	F-0227.03	Iol implantation forceps, titanium
10	H-0105	Barraquer needle holder curved, titanium
11	K-0002	Capsulotome
12	E-0000.02	Hirshman iris hook and lens manipulator, angled
13	E-0000.05	Osher & Y[hook, angled
14	T-0001	Iol spatula notched
15	T-0002	Double eye spatula
16	B-0000	Castroviejo blade holder straight, titanium
17	R-0002	Ophtalmic speculums, titanium

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant



1



2



3



4



5



6



7

8

9



10



11




12



13



14



15

16

17

OPHTHALMIC SURGICAL SETS

SURGICAL SETS OPHT

HALMIC



AO-2

FOREIGN BODY EXTRACTION

MICROSURGICAL SET

1	S-0075	Vannas scissors curved, stainless steel
2	S-0064	Castroviejo universal scissors curved, s.s.
3	S-0066	Scissors for conjunctiva curved, s.s.
4	S-0072	Troutman-Castroviejo corneal section scissors right, s.s.
5	S-0073	Troutman-Castroviejo corneal section scissors left, s.s.
6	V-013K	Inta-ocular scissors, straight
7	F-0008	Catalano tying forceps curved, titanium
8	F-0011	Catalano corneal forceps, titanium
9	F-0224	Troutman-Barraquer corneal forceps colibri type, titanium
10	F-0109	Anis anterior capsulotomy forceps, titanium
11	V-004.01T	Handle, titanium
12	V-047K	Three-prong forceps
13	V-048K	Double cup forceps
14	V-045K	Endface forceps normally close
15	F-0200	Wire removal forceps straight, titanium
16	F-0200.01	Wire removal forceps angled, titanium
17	H-0089	Tennant needle holder curved, titanium
18	B-0000	Castroviejo blade holder straight, titanium
19	E-0001	Microhook right, angled
20	E-0001.01	Microhook left, angled
21	T-0000	Eye spatula angled 0,5 mm wide
22	T-0000.01	Eye spatula angled 1,0 mm wide
23	T-0000.02	Eye spatula angled 1,5 mm wide
24	R-0001	Murdock eye speculum right, titanium
25	R-0001.01	Murdock eye speculum left, titanium
26	N-0000	Handle for cannulae, titanium
27	C-0005	Aspiration cannula 1,0 mm
28	C-0005.01	Aspiration cannula 1,2 mm
29	C-0005.02	Aspiration cannula 1,6 mm
30	C-0005.03	Aspiration cannula 2,0 mm
31		Sterilizing box

OPHTHALMIC INSTRUMENT

S



Kazan Medical Instruments Plant



1



2



3



4



5



6



7



8



9



10



11



12



13



14



15



16

OPHTHALMIC SURGICAL SETS

SURGICAL SETS OPHT

HALMIC



17



18



19



20



21



22



23



24



25



26



27



28

29

30



AO-3

IOL IMPLANTATION MICROSURGICAL SET


1	S-0075	Vannas scissors curved, stainless steel
2	S-0064	Castroviejo universal scissors curved, s.s.
3	S-0066	Scissors for conjunctiva curved, s.s.
4	S-0072	Troutman-Castroviejo corneal section scissors right, s.s.
5	S-0073	Troutman-Castroviejo corneal section scissors left, s.s.
6	S-0074	Vannas scissors angled to flat, ss
7	F-0224	Troutman-Barraquer corneal forceps colibri type, titanium
8	F-0257	Bonn iris forceps straight, titanium
9	F-0116	Terson extra capsular forceps, titanium
10	F-0109	Anis anterior capsulotomy forceps, titanium
11	F-0115	Tenagulum forceps for iol implantation, titanium
12	F-0195	Schmidt-Mumm iol holding forceps, titanium
13	F-0196	Iol implantation forceps, titanium
14	F-0197	Sergienko iol holding forceps, titanium
15	H-0089	Tennant needle holder curved, titanium
16	B-0000	Castroviejo blade holder straight, titanium
17	E-0000	Graether collar button hook, straight
18	E-0000.01	Sinskey iris and iol hook, angled
19	E-0000.02	Hirshman iris hook and lens manipulator, angled
20	E-0000.03	Bonn iris hook, angled
21	E-0000.04	Rentsch lens manipulator, straight
22	E-0000.05	Osher & Y[hook, angled
23	T-0000	Eye spatula angled 0,5 mm wide
24	T-0000.01	Eye spatula angled 1,0 mm wide
25	T-0000.02	Eye spatula angled 1,5 mm wide
26	R-0002	Eye speculum, titanium
27	K-0000	Sato knife
28	R-0007	Iridal retractor
29	N-0000	Handle for cannulae, titanium
30	C-0015	Eye knife-cannula
31	C-0016	Simcoe cortex aspiration cannula with botton port
32	C-0016.01	Simcoe cortex aspiration cannula with left port
33	C-0016.02	Simcoe cortex aspiration cannula with right port
34	C-0011	Eye aspiration cannula
35	C-0017	Park irrigation cannula
36	C-0012	Simcoe suturing needle-cannula
37	C-0013	Simcoe micro hook-cannula
38	C-0014	Anis cortex aspiration cannula
39	CT-0000	Adapter from &Luer[to &Record[, titanium
40		Sterilizing box

OPHTHALMIC INSTRUMENT

S




Kazan Medical Instruments Plant



OPHTHALMIC SURGICAL SETS

SURGICAL SETS OPHT





Kazan Medical Instruments Plant

AO-4

KIT OF INSTRUMENTS IS INTENDED FOR
SURGICAL TREATING OF CORNEAL PATHOLOGIES

1	F-0198	Hoffman-polack fixation forceps angled, titanium
2	F-0223	Pierse type iris forceps colibri, titanium
3	H-0089	Tennant needle holder curved, titanium
4	M-0002	Thornton radial marker low profile, 3 blades
5	M-0002.01	Thornton radial marker low profile, 4 blades
6	M-0002.02	Thornton radial marker low profile, 5 blades
7	K-0009	Donor cornea punch for through keratoplasty (with set of trepans K-0008 - K-0008.15)
8	Z-0001	Donor cornea holder for layer keratoplasty
9	G-0003	Sterilizing box



1



2



3



4



5



6



7



8

OPHTHALMIC INSTRUMENTS NUMERICAL INDEX

qatalogue#..... Pages	qatalogue#..... Pages	qatalogue#..... Pages	qatalogue#..... Pages
AO-1..... 160	E-0013 137	F-0231.01 44	F-0264.04..... 48
AO-2 162	F-0008..... 57	F0232 21	F-0265..... 60
AO-3 164	F-0009 50	F-0232.01 21	F-0265.01 60
AO-4 166	F-0009.01..... 50	F-0232.02 21	F-0266..... 51
B-0000 154	F-0010 59	F-0233 22	F-0267 39
B-0001..... 154	F-0011..... 48	F-0233.01 22	F-0268..... 59
B-0008 154	F-0011.01 48	F-0233.02..... 22	F-0268.01 59
C-0005 13	F-0012 57	F-0234..... 21	F-0269..... 51
C-0005.01..... 13	F-0036.04 49	F-0234.01 21	F-0270 25
C-0005.02 13	F-0037 71	F-0234.02..... 21	F-0270.01..... 25
C-0005.03 13	F0037.01..... 71	F-0235 41	F-0270.02 25
C-0011..... 13	F-0038..... 71	F-0235.01 41	F-0271..... 26
C-0011.01 13	F-0038.01 71	F-0236..... 41	F-0271.01 26
C-0012 13	F-0039..... 72	F-0236.01 41	F-0271.02..... 26
C-0013 13	F-0039.01 72	F-0237 42	F-0272 26
C-0014 14	F-0109 73	F-0237.01..... 42	F-0272.01..... 26
C-0015 14	F-0113 57	F-0238 55	F-0272.02 26
C-0016 14	F-0114 65	F-0239..... 57	F-0273 49
C-0016.01 14	F-0115 67	F-0240..... 58	F-0274 49
C-0016.02..... 14	F-0163 74	F-0241 58	F-0275 49
C-0017 15	F-0173..... 36	F-0242 61	F-0279 62
C-0018 15	F-0173.01 36	F-0243..... 29	F-0311 62
C-0019 17	F-0174..... 37	F-0243.01 29	F-0312..... 62
C-0020 17	F-0174.01 37	F-0243.02..... 29	F-0313 62
C-0021 17	F-0174.02..... 37	F-0244..... 29	F-0315.01..... 50
C-0022..... 17	F-0174.03 37	F-0244.01 29	F-0316.02 63
C-0026 15	F-0176..... 65	F-0244.02..... 29	F-0316.03 63
C-0028..... 16	F-0177 64	F-0244.03..... 29	F-0317..... 63
C-0028.01 16	F-0177.01 64	F-0245..... 30	F-0318.01..... 52
C-0028.02..... 16	F-0177.02 64	F-0245.01 30	F-0318.03 52
C-0029 16	F-0177.03..... 64	F-0245.02..... 30	F-0319.04 53
C-0032..... 16	F-0178.02..... 66	F-0246..... 32	F-0319.05 53
CT-0000 12	F-0178.03..... 66	F-0247 30	F-0319.06 53
CT-0001..... 12	F-0181 31	F-0247.01..... 30	F-0319.07 53
D-0000..... 152	F-0181.01 31	F-0247.02 30	F-0319.09 53
D-0001..... 152	F-0195 66	F-0248..... 22	F-0320.03..... 46
E-0000 136	F-0196 67	F-0248.01 22	F-0320.04..... 46
E-0000.01..... 136	F-0197..... 66	F-0249..... 23	F-0320.05..... 46
E-0000.02 136	F-0198 65	F-0249.01 23	F-0321.03 27
E-0000.03 136	F-0217..... 73	F-0250..... 23	F-0321.04 27
E-0000.04 137	F-0223 59	F-0250.01 23	F-0321.05 27
E-0000.05 137	F-0224 51	F-0251..... 45	F-0323 36
E-0000.06 136	F-0225 51	F-0252 45	F-0323.01 36
E-0000.07 136	F-0226 24	F-0253 45	F-0324..... 74
E-0000.09 136	F-0226.01 24	F-0254..... 38	F-0325 70
E-0000.12..... 136	F-0226.02 24	F-0255 38	F-0326..... 45
E-0000.13..... 137	F-0227 24	F-0256..... 39	F-0327 69
E-0001 137	F-0227.01..... 24	F-0257 38	F-0328 28
E-0001.01 137	F-0227.02 24	F-0258 38	F-0328.01 28
E-0002 138	F-0228 25	F-0259..... 39	F-0328.02..... 28
E-0002.01..... 138	F-0228.01..... 25	F-0260..... 38	F-0329..... 32
E-0003 138	F-0228.02 25	F-0261 38	F-0331 63
E-0005 138	F-0229 43	F-0262 39	F-0332 61
E-0011..... 137	F-0229.01 43	F-0263..... 39	F-0335..... 75
E-0011.01 137	F-0230..... 43	F-0264..... 48	F-0336..... 75
E-0012 138	F-0230.01 43	F-0264.01 48	F-0336.04 75
E-0012.01..... 138	F-0231..... 44	F-0264.03..... 48	F-0337 46

OPHTHALMIC INSTRUMENTS NUMERICAL INDEX

qatalogue#.....	Pages	qatalogue#..... Pages	qatalogue#..... Pages	qatalogue#.....	Pages
F-0338.....	46	F-0481.01.....	72	H-0111.01	99
F-0339.....	42	F-0484.....	55	H-0111.02	100
F-0340.....	60	F-0486.....	54	H-0112.....	98
F-0341	40	F-0495.....	35	H-0112.01	99
F-0342.....	53	F-0495.01	35	H-0112.02	100
F-0344.....	52	F-0499.....	33	H-0113	92
F-0345.....	52	F-0500.....	33	H-0113.01.....	93
F-0348.....	67	F-0500.01	33	H-0113.02	94
F-0349.....	67	F-0500.02	34	H-0114	104
F-0350.....	68	F-0501	34	H-0115	95
F-0352	68	F-0502.....	34	H-0115.01.....	95
F-0354.....	56	F-0503.....	61	H-0115.02	95
F-0355.....	56	F-0561	47	H-0115.03	95
F-0356.....	54	F-0563.....	28	H-0117	92
F-0357	54	F-0565.....	47	H-0117.01	93
F-0379	40	F-0567	50	H-0117.02.....	94
F-0384.....	60	F-0570	68	H-0119	105
F-0388.....	74	H-0055	111	H-0119.01.....	105
F-0389.02.....	27	H-0056	111	H-0121.....	98
F-0394.01	31	H-0056.01	111	H-0121.01	99
F-0395.....	31	H-0064.....	111	H-0121.02	100
F-0395.01	31	H-0064.01	111	H-0122	92
F-0395.02.....	31	H-0084	105	H-0122.01	93
F-0396.....	65	H-0089	92	H-0122.02	94
F-0397	51	H-0089.01	93	H-0122.06.....	92
F-0397.01	51	H-0089.02.....	94	H-0122.07	93
F-0422	32	H-0090.....	96	H-0122.08.....	94
F-0422.01	32	H-0090.01	96	H-0123.....	92
F-0423.....	35	H-0090.02.....	97	H-0123.01	93
F-0423.01	35	H-0100.....	96	H-0123.02.....	94
F-0428.....	42	H-0100.01	96	H-0123.06.....	92
F-0432.....	75	H-0100.02	97	H-0123.07	93
F-0434.....	56	H-0101	98	H-0123.08.....	94
F-0435.....	44	H-0101.01.....	99	H-0126.01	112
F-0435.01	44	H-0101.02	100	H-0129.01	112
F-0435.02.....	44	H-0102.....	98	H-0130.01	112
F-0436.....	55	H-0102.01	99	H-0135.....	109
F-0436.01	55	H-0102.02.....	100	H-0135.01	109
F-0437	68	H-0103.....	98	H-0136.....	108
F-0438.....	61	H-0103.01	99	H-0138.....	108
F-0457	69	H-0103.02.....	100	H-0173	101
F-0457.01	69	H-0104.....	98	H-0173.01	101
F-0458.....	70	H-0104.01	99	H-0173.02	101
F-0458.01	70	H-0104.02	100	H-0178	106
F-0458.04	70	H-0105.....	102	H-0178.01.....	106
F-0458.05	70	H-0105.01	102	H-0178.02	107
F-0459.....	40	H-0105.02.....	102	H-0178.03.....	106
F-0459.01	40	H-0106.....	103	H-0178.04.....	106
F-0459.02.....	40	H-0106.01	103	H-0178.05.....	107
F-0459.03	40	H-0106.02	103	H-0182	110
F-0460.01	47	H-0107	97	H-0183.....	110
F-0460.02	47	H-0108.....	97	H-0184.....	108
F-0461	73	H-0109.....	104	H-0212	109
F-0477	54	H-0110.....	98	H-0214.....	107
F-0477.01.....	54	H-0110.01.....	99	K-0000	141
F-0479	69	H-0110.02	100	K-0001.....	139
F-0481	72	H-0111	98	K-0002	139
				K-0003	139
				K-0005	139
				K-0006	157
				K-0007	157
				K-0008	157
				K-0009	156
				K-0010.04	139
				K-0010.05	139
				K-0010.06	139
				K-0016.....	141
				K-0017	141
				K-0018.....	141
				K-0019.....	141
				K-0021.01	8
				K-0021.02.....	8
				K-0021.03.....	8
				K-0021.04.....	8
				K-0021.05.....	8
				K-0021.06.....	8
				K-0021.07	8
				K-0021.08.....	8
				K-0025.01.....	9
				K-0025.02	9
				K-0025.03	9
				K-0025.04	9
				K-0025.05	9
				K-0025.06	9
				K-0025.07	9
				K-0025.08	9
				K-0026.01.....	7
				K-0026.02	7
				K-0026.03	7
				K-0026.04	7
				K-0026.05	7
				K-0026.06	7
				K-0026.07	7
				K-0026.08	7
				K-0026.09	7
				K-0026.11	7
				K-0026.12	7
				K-0026.15.....	7
				K-0026.21	7
				K-0027.01	6
				K-0027.02.....	6
				K-0027.03	6
				K-0027.04	6
				K-0027.05	6
				K-0027.06	6
				K-0027.07.....	6
				K-0027.08	6
				K-0027.09	6
				K-0027.10	6
				K-0027.11.....	6
				K-0027.12	6
				K-0027.13	6
				K-0031.....	140
				K-0031.01	140

OPHTHALMIC INSTRUMENTS NUMERICAL INDEX

qatalogue#..... Pages	qatalogue#..... Pages	qatalogue#..... Pages	qatalogue#..... Pages
K-0031.02..... 140	SMP-0015..... 132	S-0163.03 125	V-050 K..... 84
K-0031.03..... 140	S-0000 115	S-0164 122	V-053 K..... 81
K-0031.04 140	S-0000.01..... 115	S-0165 114	V-054 K..... 83
L-0000 142	S0018..... 124	S-0171 125	V-077 K 155
L-0002..... 142	S-0026..... 116	T-0000 147	W-0000 150
M-0001 155	S-0027 116	T-0000.01..... 147	Z-0000 149
M-0001.01 155	S-0029..... 116	T-0000.02 147	Z-0000.01..... 149
M-0002..... 142	S-0040..... 119	T-0001 147	Z-0000.02 149
M-0003 142	S-0040.01 119	T-0002..... 148	Z-0000.03 149
M-0007 143	S-0040.02 119	T-0003 148	Z-0000.07 151
M-0008 144	S-0043..... 119	T-0006 147	Z-0000.08 151
M-0009 144	S-0043.01 119	T-0006.01..... 147	Z-0000.09 151
M-0010 144	S-0043.02..... 119	T-0006.02 147	Z-0001 156
N-0000..... 12	S-0044..... 123	T-0006.03 147	Z-0005 151
N-0001..... 12	S-0044.01 123	T-0007..... 148	Z-0006 150
R-0000 128	S-0045..... 123	T-0014 146	Z-0006.01-20 .. 150
R-0001..... 128	S-0045.01 123	T-0014.01..... 146	Z-0007..... 151
R-0001.01 128	S-0046..... 123	V-004.01T..... 78	Z-0008 143
R-0002 128	S-0046.01 123	V-006 A-K..... 78	
R-0003 129	S-0047 123	V-006.01 A-K .. 78	
R-0005 133	S-0047.01 123	V-007 A-K 78	
R-0006 145	S-0064..... 117	V-007.01 A-K... 78	
R-0006.01 145	S-0065..... 124	V-008 K..... 79	
R-0006.02 145	S-0066..... 118	V-009 K..... 79	
R-0006.03..... 145	S-0067 124	V-010 K 79	
R-0006.04..... 145	S-0068..... 124	V-011 K 80	
R-0007..... 146	S-0069..... 121	V-012 K..... 80	
R-0008 133	S-0069.01 121	V-013 K 79	
R-0008.01..... 133	S-0070 124	V-014 K 79	
R-0010..... 130	S-0072 120	V-015 K 80	
R-0010.01 130	S-0073 120	V-016 K 80	
R-0010.02..... 130	S-0074 118	V-017 K..... 79	
R-0010.03..... 130	S-0075 118	V-018 K 80	
R-0010.04 130	S-0079 118	V-020 K..... 81	
R-0010.05..... 130	S-0080..... 121	V-021 K..... 81	
R-0011..... 131	S-0081 121	V-022 K 81	
R-0011.01 131	S-0122..... 118	V-023 K..... 81	
R-0012 130	S-0130 116	V-024 K..... 82	
R-0012.01 130	S-0130.01..... 116	V-025 K..... 82	
R-0013..... 131	S-0134 120	V-026 K..... 82	
R-0013.01 131	S-0134.01..... 120	V-027 A-K..... 85	
R-0013.02..... 131	S-0135 122	V-028 A-K 85	
R-0013.03..... 131	S-0135.01..... 122	V-029 A-K 85	
R-0027 133	S-0136 117	V-030 A-K 85	
R-0027.01 133	S-0137..... 117	V-031 K 84	
R-0027.02..... 133	S-0139 122	V-033 K..... 87	
R-0031..... 134	S-0140 115	V-034 K..... 87	
R-0032 132	S-0140.01..... 115	V-035 K..... 86	
R-0033 132	S-0141 117	V-036 K..... 87	
RA-100..... 134	S-0143 116	V-042 K..... 83	
RA-100.01 134	S-0146 122	V-043 K..... 83	
RT-0001..... 129	S-0146.01..... 122	V-044 K..... 83	
RT-0001.01 129	S-0160 114	V-045 K..... 84	
RT-0002 129	S-0163 125	V-047 K 86	
RT-0002.01..... 129	S-0163.01..... 125	V-048 K..... 86	
RT-0003 129	S-0163.02 125	V-049 K..... 84	

MICROSURGICAL
INSTRUMENTS CJSC
12, Saidashev str., Kazan,
Tatarstan, 420022
Russia
Tel/Fax: 007-8432-933-323
007-8432-782-805
E-mail: micro@bancorp.ru
micro@kzn.ru
www.micro-instruments.com

