

INTRA OCULAR INSTRUMENTS

INTRA OCULAR

INSTRUMENTS



Kazan Medical Instruments Plant





HANDPIECE

V-004.01T

HANDLE, TITANIUM

Overall length 102 mm; 4,0 inches


Round handle 8 mm dia



L=30mm; D=0,9mm
B=0,7mm; H=0,7mm

V-006 A-K


LONG TIP VITREOTOMY



L=18mm; D=0,9mm
B=0,7mm; H=0,7mm

V-006.01 A-K


SHORT TIP VITREOTOMY



L=30mm; D=1,2mm
B=0,9mm; H=0,9mm

V-007 A-K

LONG TIP LENSEKTOME



L=18mm; D=1,2mm
B=0,9mm; H=0,9mm


V-007.01 A-K

SHORT TIP LENSEKTOME



Kazan Medical Instruments Plant


INTRA OCULAR INSTRUMENTS



L=30mm; D=0,9mm
B=2,4mm; H=1,4mm

V-013 K


SCISSORS STRAIGHT



L=18mm; D=0,65mm
B=1,6mm; H=0,8mm

V-017 K

**ANTERIOR CAPSULOTOMY
SCISSORS STRAIGHT**



L=30mm; D=0,9mm
B=1,35mm; H=1,4mm

V-008 K

**SCISSORS VERTICAL BLUNT TIPS,
GUILLotine ACTION**



L=30mm; D=0,9mm
B=1,35mm; H=1,4mm

V-009 K


**SCISSORS VERTICAL SHARP
TIPS, GUILLotine ACTION**



L=30mm; D=0,9mm
B=1,35mm; H=1,4mm

V-010 K

**SCISSORS VERTICAL
GUILLotine ACTION,
FRONT TIP WITH OLIVE**



L=30mm; D=0,9mm
B=1,7mm; H=1,5mm



V-014 K

**SCISSORS HORIZONTAL
SHARP TIPS**

INTRA OCULAR INSTRUMENTS




Kazan Medical Instruments Plant



L=18mm; D=0,65mm
B=1,1mm; H=1,2mm

V-015 K


**ANTERIOR CAPSULOTOMY
SCISSORS, CURVED LEFT**



L=18mm; D=0,65mm
B=1,1mm; H=1,2mm

V-016 K


**ANTERIOR CAPSULOTOMY
SCISSORS, CURVED RIGHT**



L=18mm; D=0,65mm
B=1,1mm; H=1,2mm

V-018 K


**ANTERIOR CAPSULOTOMY
SCISSORS**



L=30mm; D=0,9mm
B=1,0mm; H=1,4mm

V-011 K

**SCISSORS VERTICAL STEEP
BLANT TIPS, GUILLOTINE ACTION**





L=30mm; D=0,9mm
B=1,0mm; H=1,4mm

V-012 K

**SCISSORS VERTICAL STEEP
SHARP TIPS, GUILLOTINE ACTION**




Kazan Medical Instruments Plant

INTRA OCULAR
INSTRUMENTS

L=30mm; D=0,9mm
B=1,8mm; H=1,4mm

V-020 K


**FORCEPS STRAIGHT,
SMOOTH JAWS**



L=30mm; D=0,9mm
B=1,8mm; H=1,4mm

V-021 K


**FORCEPS STRAIGHT
WITH PLATFORM,
SMOOTH JAWS**



L=30mm; D=0,9mm
B=1,8mm; H=1,4mm

V-022 K


**FORCEPS STRAIGHT,
SERRATED JAWS**



L=30mm; D=0,9mm
B=2,4mm; H=1,5mm

V-023 K

**FORCEPS STRAIGHT,
SMOOTH JAWS**



L=30mm; D=0,9mm
B=1,8mm; H=1,4mm


V-053 K

**FORCEPS STRAIGHT
WITH MICROTOOTH,
SMOOTH JAWS**

INTRA OCULAR INSTRUMENTS




Kazan Medical Instruments Plant



L=30mm; D=0,9mm
B=1,2mm; H=1,3mm


V-024 K



L=30mm; D=0,9mm
B=1,2mm; H=1,3mm

V-024 K

**FORCEPS CURVED,
SMOOTH JAWS**



L=30mm; D=0,9mm
B=1,2mm; H=1,3mm

V-025 K

**FORCEPS CURVED,
SERRATED JAWS**



L=30mm; D=0,9mm
B=1,7mm; H=1,3mm



V-026 K

**FORCEPS CURVED,
SMOOTH JAWS**



Kazan Medical Instruments Plant


INTRA OCULAR INSTRUMENTS



L=30mm; D=1,0mm
B=1,0mm; H=1,0mm

V-042 K


**FORCEPS SPHERICAL
JAWS,
NORMALLY OPEN**



L=30mm; D=1,0mm
B=1,0mm; H=1,0mm

V-043 K

**FORCEPS SPHERICAL
JAWS,
NORMALLY CLOSE**



L=25mm; D=1,2mm
B=2,7mm; H=1,0mm

V-044 K

**MICROSURGICAL
NEEDLE HOLDER**

● TC jaws



L=30mm; D=0,9mm
B=1,6mm; H=1,4mm


V-054 K

CAPSULORHEXIS FORCEPS

INTRA OCULAR INSTRUMENTS




Kazan Medical Instruments Plant



L=30mm; D=0,9mm
B=2,2mm; H=2,0mm

V-031 K


**FORCEPS SIDE GRIPPING
STEEP SMOOTH JAWS**



L=30mm; D=1,2mm
B=1,2mm; H=2,5mm

V-045 K

**ENDFACE FORCEPS
NORMALLY CLOSE**



L=30mm; D=0,9mm
B=1,65mm; H=3,0mm

V-049 K

FORCEPS ANGLED SMALL

- Normally close
- Serrated jaws



L=30mm; D=1,2mm
B=2,4mm; H=3,5mm


V-050 K

FORCEPS ANGLED LARGE

- Normally close
- Serrated jaws




Kazan Medical Instruments Plant



L=30mm; D=1,0mm
B=1,2mm; H=1,5mm

V-027 A-K

**ASPIRATION FORCEPS
SHARP SMALL JAWS**



L=30mm; D=1,0mm
B=2,0mm; H=1,5mm

V-028 A-K

**ASPIRATION FORCEPS
SHARP LARGE JAWS**



L=30mm; D=1,0mm
B=1,2mm; H=1,5mm

V-029 A-K

**ASPIRATION FORCEPS
BLUNT SMALL JAWS**



L=30mm; D=1,0mm
B=2,0mm; H=1,5mm



V-030 A-K

**ASPIRATION FORCEPS
BLUNT LARGE JAWS**

INTRA OCULAR INSTRUMENTS



Kazan Medical Instruments Plant




L=30mm; D=1,2mm
B=1,4mm; H=4,0mm

V-048 K

DOUBLE CUP FORCEPS

● Normally close




L=30mm; D=1,2mm
B=3,5mm; H=4,0mm

V-035 K

TWO-LEAF FORCEPS

● Normally close



L=30mm; D=1,2mm
B=4,8mm; H=4,5mm



V-047 K

THREE-PRONG FORCEPS

● Normally close




Kazan Medical Instruments Plant

INTRA OCULAR
INSTRUMENTS

$L=18\text{mm}$; $D=0,65\text{mm}$
 $B=1,2\text{mm}$; $H=0,8\text{mm}$

V-033 K


ANTERIOR CAPSULE
 FORCEPS CURVED



$L=18\text{mm}$; $D=0,65\text{mm}$
 $B=1,35\text{mm}$; $H=2,5\text{mm}$

V-036 K

DOUBLE HOOK FORCEPS
 CURVED



$L=18\text{mm}$; $D=0,65\text{mm}$
 $B=0,9\text{mm}$; $H=1,5\text{mm}$

V-034 A-K

CAPSULORHEXIS
 FORCEPS WITH
 ASPIRATION

● One sharp jaw

OPHTHALMIC INSTRUMENTS



Kazan Medical Instruments Plant

**NEEDLE
HOLDERS**

TITANIUM




UNIVERSAL LOCK FOR NEEDLE HOLDERS

NEW

*U[®] Universal needle holder can be used as a needle holder both **with** the lock and **without** it. The removable lock allows a surgeon the option to use or not to use the lock during operation.*

U[®] Универсальный иглодержатель может быть использован в качестве иглодержателя с замком, так и без него. Съемный замок дает хирургу возможность использовать замок или не использовать во время операции.


Figure 1
Position **with lock**
Позиция **с замком**




 **A - to move the catch to **without** lock position**
А - переместить защелку в положение без замка

UNIVERSAL LOCK FOR NEEDLE HOLDERS

Figure 2
Position without lock
Позиция без замка



 **B - to move the pin to with lock position**
В - переместить упор в положение с замком



TC - Tungsten Carbide

Tungsten Carbide particles metallurgically fused onto the jaws create a fine textured surface. This extremely durable surface allows for the secure, atraumatic handling of delicate sutures and tissues.

ТС - Карбид вольфрама

Карбид вольфрама, напыленный на губки инструмента, создает поверхность с мелкой текстурой.

Это чрезвычайно прочная поверхность обеспечивает надежную атравматичную работу с деликатными шовным материалом и тканями.

NEEDLE HOLDERS



**NEEDLE HOLDER,
TITANIUM**

- Overall length 135 mm; 5,3 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11 mm long

H-0089

TENNANT



Tips 0,6 mm wide
Jaws 7 mm long

H-0113

BARRAQUER



Tips 0,4 mm wide
Jaws 11 mm long

H-0117

CASTROVIEJO



Tips 0,6 mm wide
Jaws 10 mm long

H-0122

CASTROVIEJO



Tips 0,6 mm wide
Jaws 10 mm long

H-0122.06

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0123

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0123.06

CASTROVIEJO



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 135 mm; 5,3 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11 mm long

H-0089.01

TENNANT



Tips 0,6 mm wide
Jaws 7 mm long

H-0113.01

BARRAQUER



Tips 0,4 mm wide
Jaws 11 mm long

H-0117.01

CASTROVIEJO



Tips 0,6 mm wide
Jaws 10 mm long

H-0122.01

CASTROVIEJO



Tips 0,6 mm wide
Jaws 10 mm long

H-0122.07

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0123.01

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0123.07

CASTROVIEJO

NEEDLE HOLDERS



**NEEDLE HOLDER
WITH UNIVERSAL LOCK,
TITANIUM**

- Overall length 135 mm; 5,3 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11 mm long

H-0089.02

TENNANT



Tips 0,6 mm wide
Jaws 7 mm long

H-0113.02

BARRAQUER



Tips 0,4 mm wide
Jaws 11 mm long

H-0117.02

CASTROVIEJO



Tips 0,6 mm wide
Jaws 10 mm long

H-0122.02

CASTROVIEJO



Tips 0,6 mm wide
Jaws 10 mm long

H-0122.08

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0123.02

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0123.08

CASTROVIEJO



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



Tips 0,4 mm wide
Jaws 11 mm long

H-0115

STRAIGHT



Tips 0,4 mm wide
Jaws 11 mm long

H-0115.02

CURVED

McPHERSON NEEDLE HOLDER, TITANIUM

- Overall length 110 mm; 4,3 inches
- Round handle 6 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11 mm long

H-0115.01

STRAIGHT



Tips 0,4 mm wide
Jaws 11 mm long

H-0115.03

CURVED

McPHERSON NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 110 mm; 4,3 inches
- Round handle 6 mm dia
- TC jaws

NEEDLE HOLDERS



Tips 0,4 mm wide
Jaws 7,5 mm long

H-0090
COHAN



Tips 0,4 mm wide
Jaws 8,5 mm long

H-0100
BARRAQUER

NEEDLE HOLDER, TITANIUM

- Overall length 100 mm; 3,9 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 7,5 mm long

H-0090.01
COHAN



Tips 0,4 mm wide
Jaws 8,5 mm long

H-0100.01
BARRAQUER


NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 100 mm; 3,9 inches
- Round handle 8 mm dia
- TC jaws



Kazan Medical Instruments Plant


NEEDLE HOLDERS NEEDLE HOLDERS



Tips 0,4 x 0,4 mm wide
Jaws 7,5 mm long

H-0090.02

COHAN



Tips 0,4 x 0,35 mm wide
Jaws 8,5 mm long

H-0100.02

BARRAQUER

NEEDLE HOLDER WITH UNIVERSAL LOCK, TITANIUM


- Overall length 100 mm; 3,9 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 10 mm long

H-0107

BARRAQUER



Tips 0,4 mm wide
Jaws 8 mm long

H-0108

BECHERT-SINSKEY

NEEDLE HOLDER, TITANIUM


- Overall length 100 mm; 3,9 inches
- Round handle 6 mm dia
- TC jaws

NEEDLE HOLDERS




NEEDLE HOLDER, TITANIUM


- Overall length 125 mm; 4,9 inches
- Round handle 8 mm dia
- TC jaws




CASTROVIEJO




CASTROVIEJO




BARRAQUER



CASTROVIEJO




BECHERT



McPHERSON



BARRAQUER



CASTROVIEJO



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 125 mm; 4,9 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11 mm long

H-0101.01

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0102.01

CASTROVIEJO



Tips 0,4 mm wide
Jaws 7 mm long

H-0103.01

BARRAQUER



Tips 0,4 mm wide
Jaws 7 mm long

H-0104.01

CASTROVIEJO



Tips 0,4 mm wide
Jaws 8 mm long

H-0110.01

BECHERT



Tips 0,4 mm wide
Jaws 10,5 mm long

H-0111.01

McPHERSON



Tips 0,4 mm wide
Jaws 10 mm long

H-0112.01

BARRAQUER



Tips 0,5 mm wide
Jaws 10 mm long

H-0121.01

CASTROVIEJO

NEEDLE HOLDERS



NEEDLE HOLDER WITH UNIVERSAL LOCK, TITANIUM

- Overall length 125 mm; 4,9 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11 mm long

H-0101.02

CASTROVIEJO



Tips 0,4 mm wide
Jaws 11 mm long

H-0102.02

CASTROVIEJO



Tips 0,4 mm wide
Jaws 7 mm long

H-0103.02

BARRAQUER



Tips 0,4 mm wide
Jaws 7 mm long

H-0104.02

CASTROVIEJO



Tips 0,4 mm wide
Jaws 8 mm long

H-0110.02

BECHERT



Tips 0,4 mm wide
Jaws 10,5 mm long

H-0111.02

McPHERSON



Tips 0,4 mm wide
Jaws 10 mm long

H-0112.02

BARRAQUER



Tips 0,5 mm wide
Jaws 10 mm long

H-0121.02

CASTROVIEJO



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



Tips 0,4 mm wide
Jaws 11,5 mm long

H-0173

Overall length 135 mm; 5,3 inches



Tips 0,55 mm wide
Jaws 11,5 mm long

H-0173.02

Overall length 145 mm; 5,7 inches

NEEDLE HOLDER, TITANIUM

- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 11,5 mm long

H-0173.01



NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 135 mm; 5,3 inches
- Round handle 8 mm dia
- TC jaws

● All needle holders also available in stainless steel

● TC-Tungsten Carbide

NEEDLE HOLDERS



Tips 0,4 mm wide
Jaws 10 mm long



H-0105



H-0105.01
with lock



H-0105.02
with universal lock

BARRAQUER NEEDLE HOLDER, TITANIUM

- Overall length 125 mm; 4,9 inches
- Round handle 7 mm dia
- TC jaws

● All needle holders also available in stainless steel

TC-Tungsten Carbide




Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



H-0106



H-0106.01
with lock



H-0106.02
with universal lock

TENNAIT NEEDLE HOLDER CURVED, TITANIUM

- Overall length 135 mm; 5,3 inches
- Round handle 10 mm dia
- TC jaws

NEEDLE HOLDERS



Tips 0,4 mm wide
Jaws 11 mm long

H-0114



BARRAQUER NEEDLE HOLDER CURVED, TITANIUM

- Overall length 135 mm; 5,3 inches
- Round handle 6 mm dia
- TC jaws



Tips 0,4
Jaws 8 mm long

H-0109



BECHERT-SINSKEY NEEDLE HOLDER STRAIGHT, TITANIUM

- Overall length 125 mm; 4,9 inches
- Round handle 6 mm dia
- TC jaws



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



Tips 0,4 mm wide
Jaws 8,5 mm long

H-0119

STRAIGHT



Tips 0,4 mm wide
Jaws 8,5 mm long

H-0119.01

CURVED

NEEDLE HOLDER, TITANIUM

- Overall length 105 mm; 4,1 inches
- Round handle 7 mm dia with two hinge joints
- TC jaws



Tips 0,8 mm wide
Jaws 10 mm long

H-0084

KALT NEEDLE HOLDER STRAIGHT, TITANIUM

- Overall length 135 mm; 5,3 inches
- Flat locking handle 10 mm wide
- TC jaws

● All needle holders also available in stainless steel

● TC-Tungsten Carbide

NEEDLE HOLDERS



Tips 0,8 wide
Jaws 9,5 mm long

H-0178



Tips 0,8 mm wide
Jaws 9,0 mm long

H-0178.03

NEEDLE HOLDER, TITANIUM

- Overall length 140 mm; 5,5 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,8 mm wide
Jaws 9,5 mm long

H-0178.01



Tips 0,8 mm wide
Jaws 9,0 mm long

H-0178.04

NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 140 mm; 5,5 inches
- Round handle 8 mm dia
- TC jaws



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



Tips 0,8 mm wide
Jaws 9,5 mm long

H-0178.02



Tips 0,8 mm wide
Jaws 9,0 mm long

H-0178.05

NEEDLE HOLDER WITH UNIVERSAL LOCK, TITANIUM

- Overall length 140 mm; 5,5 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,9 mm wide
Jaws 11 mm long

H-0214

NEEDLE HOLDER, TITANIUM

- Overall length 135 mm; 5,3 inches
- Round handle 8 mm dia
- TC jaws

NEEDLE HOLDERS



Tips 1,5 mm wide
Jaws 10 mm long

H-0184



Tips 0,25 mm wide
Jaws 10 mm long

H-0138

NEEDLE HOLDER, TITANIUM

- Overall length 135 mm; 5,3 inches
- Round handle 8,0 mm dia
- TC jaws



Tips 0,5 mm wide
Jaws 10 mm long

H-0136



NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 130 mm; 5,1 inches
- Round handle 8 mm dia
- TC jaws



Kazan Medical Instruments Plant

NEEDLE HOLDERS NEEDLE HOLDERS



Tips 0,6 mm wide
Jaws 11 mm long

H-0135



Tips 0,6 mm wide
Jaws 11 mm long

H-0135.01

NEEDLE HOLDER WITH LOCK, TITANIUM

- Overall length 140 mm; 5,5 inches
- Flat handle 8,5 mm wide
- TC jaws



Tips 0,5 mm wide
Jaws 8 mm long

H-0212



NEEDLE HOLDER, TITANIUM

- Overall length 110 mm; 4,3 inches
- Round handle 8,5 mm dia
- TC jaws

● All needle holders also available in stainless steel

● TC-Tungsten Carbide

NEEDLE HOLDERS



Tips 0,5 mm wide
Jaws 12 mm long

H-0182

NEEDLE HOLDER, STAINLESS STEEL

- Overall length 135 mm; 5,3 inches
- Round handle 10 mm dia
- TC jaws



Tips 1,5 mm wide
Jaws 12 mm long

H-0183

NEEDLE HOLDER, STAINLESS STEEL

- Overall length 135 mm; 5,3 inches
- Round handle 10 mm dia
- TC jaws



Kazan Medical Instruments Plant

NEEDLE HOLDER-SCISSORS



Tips 0,5 mm wide
Jaws 3,5 mm long
Blades 7,5 mm long

H-0055



Tips 0,4 mm wide
Jaws 2,5 mm long
Blades 6,5 mm long

H-0056



Tips 0,4 mm wide
Jaws 2,5 mm long
Blades 6,5 mm long

H-0064

SNEAD NEEDLE HOLDER- SCISSORS, STAINLESS STEEL

- Overall length 140 mm; 5,5 inches
- Round handle 8 mm dia
- TC jaws



Tips 0,4 mm wide
Jaws 2,5 mm long
Blades 6,5 mm long

H-0056.01



Tips 0,4 mm wide
Jaws 2,5 mm long
Blades 6,5 mm long

H-0064.01

SNEAD NEEDLE HOLDER-SCISSORS WITH LOCK, STAINLESS STEEL

- Overall length 140 mm; 5,5 inches
- Round handle 10 mm dia
- TC jaws

NEEDLE HOLDER-SCISSORS OPHTHALMIC INSTRUMENTS

HOLDER-SCISSORS




EDLE



Tips 0,5 mm wide
Jaws 3,5 mm long
Blades 7,5 mm long
H-0126.01



Tips 0,5 mm wide
Jaws 3,5 mm long
Blades 7,5 mm long
H-0130.01



Tips 0,5 mm wide
Jaws 3,5 mm long
Blades 7,5 mm long
H-0129.01

- SNEAD NEEDLE
- TC jaws
- HOLDER-SCISSORS,
- STAINLESS STEEL

Overall length 140 mm; 5,5 inches
Flat handle 8 mm wide

- SNEAD NEEDLE
- HOLDER-SCISSORS,
- STAINLESS STEEL

Overall length 105 mm; 4,1 inches
Flat handle 8 mm wide
TC jaws

NEEDLE

- HOLDER-SCISSORS,
- STAINLESS STEEL

Overall length 105 mm; 4,1 inches
Round handle 8 mm dia
TC jaws